

## Data sheet

### TIDAS S 700 / CCD UV/NIR 2098



#### Specifications:

Wavelength range	200 - 980 nm
Spectral resolution	< 4 nm
Wavelength accuracy	< 0,5 nm
Wavelength reproducibility	< 0,1 nm
Baseline drift @ 250 nm	$5 \cdot 10^{-4}$ AU/h according to ASTM E685
Signal-to-noise ratio	< $4 \cdot 10^{-5}$ AU according to ASTM E685
Number of diodes	1044 x 64
Dimensions	345 x 145 x 315 mm (W x H x D)
Weight	8 kg
Voltage	85 - 265 VAC / 47 - 63 Hz
IN/OUT	Standard: 2 x IN / 2 x OUT
Illumination	2 white LED, SMA 905
Optical connection	SMA 905
A/D converter	16 Bit
Order number	81 660 70

#### Optional with integrated light source:

With integrated D2- & halogen light source	81 660 06
--	-----------

Light up your vision