

System Requirements UVS 3.x / UVS X 3.x

Hardware Requirements

UVS 3.x

- 4 GB RAM (recommended 8 GB)
- QuadCore Prozessor (e.g. Intel® Core™ i5)
- 256 GB HDD (SSD recommended)
- min. 1x RS-232 port
- min. 1 USB Port for UVS-Dongle

UVS X 3.x – 21 CFR Part 11

- 8 GB RAM
- QuadCore Prozessor (e.g. Intel® Core™ i5)
- 256 GB HDD (SSD recommended)
- min. 1x RS-232 port
- min. 1 USB Port for UVS Dongle
- Gigabit-LAN

Software Requirements / IT-Environment

UVS 3.x

- Operating System Windows 10 (Professional or Enterprise Edition)
- available Windows-User with (temporary) administrator rights
- local user (if domain login is required, domain login must comply with UVS X 3 - 21 CFR Part 11 requirements)

UVS X 3.x – 21 CFR Part 11

- **currently:** Operation System Windows 10 Enterprise LTSB 2016
starting October 2021: Windows 10 Enterprise LTSB 2016 and LTSC 2019
- available Windows-User with (temporary) administrator rights
- SQL Database-Server MS SQL 2016 or 2019, Standard or Enterprise Edition – other database servers on request
A local installation will be performed using MS SQL Express.
- SQL-user with (temporary) administrator rights)
- Local user-login and domain-based user-login are supported